# **DEUTERIUM ARC LAMP - P721**

# UV Light Source for Spectroscopy - 3 volt

#### Instrument

#### **Agilent (Hewlett Packard)**

#### **Electrical**


Filament/Heater Voltage Start 3.0 VDC
Operating 1.5~2.0 VDC
Filament/Heater Current Warm-up 4.5A Max
Discharge/Arc Voltage Strike @25° > 350VDC

Operating 60 - 90 VDC

Discharge/Arc Current 10 - 300mA

Leads Black Filament & Cathode Red Anode

Terminals 3-Pin Housing



Dimensions in millimetres

#### **Optical**

Window UV Transmitting Borosilicate

Aperture type 1mm Circular Wavelength range 185nm - 400nm

Noise  $< 1 \times 10^{-4}$  Absorption units Drift  $< 3 \times 10^{-3}$  Absorption units

# Operational life

When radiant energy reaches 50% of the initial energy. Typical life 1500 hours


This Deuterium lamps is guaranteed for a period of one year from date of delivery or 1000 hours usage, which ever comes first.

### Warranty

Free of charge replacement lamp will be provided where Photron is satisfied the following conditions have been adhered to.

- Lamp is returned to supplier freight paid.
- · Lamp has been operated within our recommended specifications.
- Lamp has not failed due to accidental, misuse or transit damage.
- If a replacement lamp is supplied the validity of guarantee shall commence from delivery date of original lamp.

# Spectrum of Deuterium Emissions



#### Precautions for use

- Do not subject the lamp to mechanical shock.
- Ensure glass envelope is free from foreign material. Lamp may be wiped using alcohol liquid.
- · Caution! UV transmitting window is very thin. Handle with care.
- Care must be taken not to contact lamp with bare hand, temperature can reach over 200°C during operation.
- UV Radiation is harmful to the eyes. Do not observe light output directly.
- Deuterium lamps operate at high voltages. Extreme caution is required to prevent electrical shock.



Unit 5 - 3 Vesper Drive, Narre Warren, Victoria, Australia 3805

Tel: 6-13-9704-9933 • Fax: 6-13-9704-6289

sales@photron.com.au • service@photron.com.au

www.photronlamp.com